RESEARCH Open Access

Two target genes based multiple cross displacement amplification combined with a lateral flow biosensor for the detection of *Mycobacterium tuberculosis* complex

Junfei Huang¹, Ziyu Xiao^{1,2}, Xinggui Yang^{1,2}, Xu Chen³, Xiaojuan Wang^{1,2}, Yijiang Chen¹, Wenlin Zheng¹, Wei Chen¹, Huijuan Chen¹ and Shijun Li^{1,2*}

Abstract

Background: Tuberculosis (TB) is a serious chronic infectious disease caused by *Mycobacterium tuberculosis* complex (MTBC). Hence, the development of a novel, simple, rapid and sensitive method to detect MTBC is of great significance for the prevention and treatment of TB.

Results: In this study, multiple cross displacement amplification (MCDA) combined with a nanoparticle-based lateral flow biosensor (LFB) was developed to simultaneously detect two target genes (*IS6110* and *mpb64*) of MTBC (MCDA-LFB). One suite of specific MCDA primers designed for the *IS6110* and *mpb64* genes was validated using genomic DNA extracted from the reference strain H37Rv. The MCDA amplicons were analyzed using a real-time turbidimeter, colorimetric indicator (malachite green, MG) and LFBs. The optimal amplification temperature and time were confirmed, and the MCDA-LFB method established in the current report was evaluated by detecting various pathogens (i.e., reference strains, isolates and clinical sputum samples). The results showed that the two sets of MCDA primers targeting the *IS6110* and *mpb64* genes could effectively detect MTBC strains. The optimal reaction conditions for the MCDA assay were determined to be 67 °C for 35 min. The MCDA assay limit of detection (LoD) was 100 fg per reaction for pure genomic DNA. The specificity of the MCDA-LFB assay was 100%, and there were no cross-reactions for non-MTBC strains. For sputum samples and MTBC strain detection, the positive rate of MCDA-LFB for the detection of MTBC strains was consistent with seminested automatic real-time PCR (Xpert MTB/RIF) and higher than acid-fast staining (AFS) and culture assays when used for sputum samples. The MCDA-LFB assay was a rapid tool, and the whole procedure for MCDA-LFB, including DNA template preparation, MCDA reaction and amplification product analysis, was completed within 70 min.

Conclusion: The MCDA-LFB assay targeting the *IS6110* and *mpb64* genes is a simple, rapid, sensitive and reliable detection method, and it has potential significance for the prevention and treatment of TB.

Keywords: *Mycobacterium tuberculosis* complex, Multiple cross displacement amplification, Lateral flow biosensor, Detection, Diagnosis

*Correspondence: zjumedjun@163.com

Background

Tuberculosis (TB) is a chronic infectious disease caused by *Mycobacterium tuberculosis* complex (MTBC), and the MTBC pathogens evolve *M. tuberculosis*, *M. bovis*, *M.*



© The Author(s) 2021. **Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and the use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit http://creativeccommons.org/licenses/by/4.0/. The Creative Commons Public Domain Dedication waiver (http://creativecommons.org/publicdomain/zero/1.0/) applies to the data made available in this article, unless otherwise stated in a credit line to the data.

² Public Health School, Guizhou Medical University, Guiyang 550025, Guizhou, People's Republic of China Full list of author information is available at the end of the article

Huang et al. BMC Microbiol (2021) 21:267 Page 2 of 11

bovis Bacillus Calmette-Guerin (BCG), M. africanum, M. carinii, M. suricattae, M. orygis, M. microti, M. caprae, M. mungi, M. canettii, M. pinnipedii, and M. vole. In particular, M. tuberculosis, M. bovis and M. africanum are highly pathogenic bacteria [1-4]. TB seriously endangers human health and is a public health and social problem of global concern [5]. The World Health Organization (WHO) has listed TB as one of the major infectious diseases. In 2018, there were approximately 10 million new TB cases and 1.5 million TB-related deaths worldwide [5]. Approximately one-third of the world's people have been and/or are now infected with MTB, according to estimates from the World Health Organization [6, 7]. Thus, quick and accurate strategies for the detection and identification of MTBC strains are important for the prevention and treatment of TB. In general, the diagnosis methods for TB usually rely on conventional sputum smear microscopy (SSM) and culture identification of the organism (namely, mycobacterial cultivation identification and biochemical tests). However, deficiencies for the abovementioned conventional examination methods (including low sensitivity, time-consuming, and complicated operation steps) cannot meet the requirements for rapid and specific detection of MTBC. Hence, simple, fast, accurate and reliable detection methods are required for the detection of MTBC in laboratory diagnosis [6, 8]. With the development of molecular techniques, polymerase chain reaction (PCR) and PCR-based assays (e.g., real-time PCR, nested PCR, and GeneXpert) have been widely used for the diagnosis of TB [9-12]. Although PCR and PCR-based assays have excellent sensitivity and reliability, the needs for specific instruments and/or reagents hinder their application in basic laboratories. To overcome the shortcomings of PCR techniques, isothermal amplification techniques, including loop-mediated isothermal amplification (LAMP), cross-priming isothermal amplification (CPA), and multiple cross displacement amplification (MCDA), were developed and applied for the diagnosis of TB [13–15]. In particular, the MCDA assay based on the isothermal strand-displacement polymerization reaction, which is a highly specific and sensitive detection technique, has been established and implemented in previous studies [16, 17]. The target sequence has ten special primers spanning ten distinct regions and requires a constant temperature to react [16–18]. Real-time turbidity, agarose gel electrophoresis, colorimetric indicators and nanoparticle-based lateral flow biosensors (LFBs) were selected for MCDA product analysis, especially MCDA combined with nanoparticlebased LFB (MCDA-LFB), which makes product analysis simple and visual [15–19]. The MCDA technique was applied for the detection of Listeria monocytogenes [18], Salmonella spp. strains and Shigella spp. [19] and Vibrio

parahaemolyticus [20]. Our team also created the MCDA method for *Brucella* spp. [17] and *Neisseria meningitidis* [16] detection successfully.

In this study, the MCDA-LFB method for MTBC detection was successfully established by our group. Two specific target genes, *IS6110* and *mpb64*, were chosen for MTBC-MCDA detection, and 10 specific primers were designed for the targets. The MCDA products were analyzed using LFB, real-time turbidity and colorimetric indicators (malachite green, MG). Then, the MCDA reaction conditions, including amplification temperature and time, were optimized. Subsequently, the sensitivity and specificity of the MTBC-MCDA-LFB technique were tested, and then MCDA-LFB was applied to detect sputum specimens and MTBC strains from clinical samples.

Results

Confirmation and detection of MCDA products

The disposable lateral flow biosensor (LFB) consists of test Line 1 (TL1), test Line 2 (TL2) and a control line (CL). First, the *IS6110* gene (Fig. 1A) and *mpb64* gene (Fig. 1B) were examined by MCDA amplification and detected with LFB and MG, respectively. Then, both target genes were detected at the same time (Fig. 1C). Genomic templates of the standard strain *M. tuberculosis* (H37Rv, ATCC 27294) were used for the MCDA assay.

Optimization of the temperatures for the MTBC-MCDA assay

To evaluate the optimum amplification temperature, *M. tuberculosis* (H37Rv, ATCC 27294) strain genomic templates were used as the positive controls at a level of 100 pg per reaction, and the reactions were monitored by the real-time turbidity (LA-320C) method. Both the *IS6110* gene and *mpb64* gene were detected, and the effect was examined at fixed temperatures ranging from 63 to 70°C with 1°C intervals for MCDA amplification. The *M. avium* genomic templates were used as negative controls. According to Table 1 and Fig. 2, at 67°C, both *IS6110* and *mpb64* gene amplification had the shortest time with higher turbidity and faster reaction times at 30 and 23 min. Thus, an amplification temperature of 67°C was applied to perform the remaining experiments in the study.

Optimization of reaction time for the MTBC-MCDA-LFB assay

The four reaction times (20, 30, 40 and 50 min) were tested at 67 °C according to the standard MCDA conditions for the optimum time by the MTBC-MCDA-LFB assay during the reaction stage. The DNA level with 100 pg of *M. tuberculosis* genomic templates per reaction was displayed by three visible-red lines (TL1,

Huang et al. BMC Microbiol (2021) 21:267 Page 3 of 11




Fig. 1 Confirmation and detection of MTBC-MCDA products. **A** The LFB and MG methods were applied for *IS6110* gene amplification. The products of the MTBC-MCDA assay were visually analyzed by observation of TL1 and color change. Tube 1/biosensor 1: positive amplification of *M. tuberculosis*; Tube 2/biosensor 2: black control of DW. **B** The LFB and MG methods were applied for *mpb64* gene amplification. The products of the MTBC-MCDA assay were visually analyzed by observation of TL2 and color change. Tube 3/biosensor 3: positive amplification of *M. tuberculosis*; Tube 4/biosensor 4: black control of DW. **C** The LFB and MG methods were applied for both *IS6110* and *mpb64* gene amplification. The products of the MTBC-MCDA assay were visually analyzed by observation of TL1, TL2 and color change. Tube 5/biosensor 5: positive amplification of *M. tuberculosis*; Tube 6/biosensor 6: black control of DW

Table 1 Reaction temperature optimization for MTBC-MCDA primers

Temperature (°C)	<i>IS6110</i> gene		mpb64 gene	
	Peak of turbidity	Time (min)	Peak of turbidity	Time (min)
63	0.16	34	0.18	29
64	0.14	31	0.20	27
65	0.17	32	0.21	26
66	0.18	30	0.21	25
67	0.18	30	0.22	23
68	0.16	34	0.21	25
69	0.15	38	0.20	25
70	0.13	55	0.20	26

The threshold of turbidity> 0.1 was judged as positive for the MCDA reaction by a real-time turbidimeter

TL2 and CL) on the LFB. The earliest test lines were observed for both TL1 and TL2 when the amplification lasted for 30 min (the threshold of turbidity>0.1) (Fig. 3). To ensure adequate amplification, a reaction time of 35 min was recommended as a reasonable reaction time for the MTBC-MCDA-LFB assay in this research.

Sensitivity of the MTBC-MCDA-LFB assay

The genomic DNA templates of M. tuberculosis (H37Rv, ATCC 27294) were serially diluted (100 ng, 10 ng, 1 ng, 100 pg, 10 pg, 1 pg, 100 fg, 10 fg and 1 fg per microliter) for MTBC-MCDA sensitivity analysis. The limiting dilution of M. tuberculosis genomic DNA was evaluated by MCDA detection. The single-MCDA-LFB assay limit of detection (LoD) targeting the IS6110 gene (Fig. 4A) or mpb64 gene (Fig. 4B) was found to be 10 fg per reaction. Double target genes were found to be 100fg per reaction (Fig. 4C). As expected, CL, TL1 and/or TL 2 were observed on the biosensor, displaying positive MCDA results for the IS6110 gene and/or mpb64 gene. Double distilled water (DW) was the template for the blank control. Moreover, the analytical sensitivity of MTBC-MCDA with the biosensor was consistent with colorimetric indicator analysis.

Specificity of the MTBC-MCDA-LFB assay

Thirty-eight MTBC strains, 12 NTM strains and 10 other bacteria were used to determine the specificity of the MCDA-LFB assay. When the genomic DNA of the bacteria listed in Table 2 was used for MCDA-LFB detection, only the DNA of the MTBC strains tested positive. The genomic DNA of NTM and other bacterial strains tested negative by the MCDA-LFB assay.

Huang et al. BMC Microbiol (2021) 21:267 Page 4 of 11




Fig. 2 Optimization of reaction temperature for MTBC-MCDA primers. Both the *IS6110* gene (**A1-A8**) and *mpb64* gene (**B1-B8**) were detected at different temperatures. The genes were detected for different temperatures (63–70 °C, 1 °C intervals) with template DNA at the level of 100 pg per reaction. The standard MCDA reactions for the detection of *M. tuberculosis* were monitored by a real-time turbidimeter, the threshold value was 0.1, and a turbidity > 0.1 was set as positive. The *M. avium* genomic templates were used as a negative control




Fig. 3 The optimal duration of time required for the MTBC-MCDA-LFB method. Four distinct reaction times (20, 30, 40, and 50 min) were examined and compared at 67 °C. Biosensors represent DNA levels of 100 pg of *M. tuberculosis* templates

Application of the MTBC-MCDA-LFB assay for sputum samples

Fifty-one sputum samples (provided by pulmonary hospital of Guiyang) were detected by acid-fast staining (AFS), conventional culture method, seminested automatic real-time PCR (Xpert MTB/RIF) and MCDA-LFB. The AFS results consisted of 26 positive samples and 25 negative samples, and the positive detection rate was 50.98%. Additionally, the sputum samples were cultured, and the results showed that 35 samples were positive and 16 were negative after 8 weeks. The Xpert MTB/RIF results indicated that 39 were positive and 12 were negative, and the rate of positive detection was 76.47%. The MCDA-LFB results showed that 41 were positive and 10 were negative, and the rate of positive detection was 80.39% (Table 3).

Discussion

To date, TB remains a chronic infectious disease caused by MTBC members that seriously endangers human health [1, 5]. Thus, convenient, rapid, sensitive and specific detection of MTBC is important for the prevention Huang et al. BMC Microbiol (2021) 21:267 Page 5 of 11

and treatment of TB, especially the development of novel rapid detection techniques. However, conventional detection methods (i.e., SSM, mycobacterial cultivation identification and biochemical tests) usually cannot meet the needs of rapid detection in terms of the detection period and sensitivity [6, 8]. PCR and PCR-based assay techniques displayed certain reliability and sensitivity in previous publications, but the requirements for PCR thermal cyclers and expensive reagents hindered their application in resource-poor areas [8]. Among the various rapid detection methods, MCDA, as a low-cost (the cost of a single MCDA reaction was estimated to be approximately 5 USD, and the LFB was 2.5 USD per test), highly sensitive and specific assay based on the LAMP technique [21], was established and applied for the detection of various pathogens (e.g., bacteria, viruses, and fungi) [21–23]. Currently, verification methods for MCDA amplicons have always been a major concern, especially for multiplex MCDA assays. Unfortunately, conventional methods (containing visual reagents, agarose gel electrophoresis and real-time turbidimetry) have difficulty verifying multiplex MCDA amplicons. Thus, LFBs based on nanoparticles were designed and utilized in our experiments. In this report, the MCDA technique combined with a nanoparticle-based LFB (MCDA-LFB) for rapid detection of MTBC was developed successfully.

In the current report, the IS6110 and mpb64 genes were chosen as specific molecular targets, and two sets of MCDA primers were designed. The IS6110 sequence, which belongs to a family of ISs of the IS3 category with high specificity, has been widely applied for MTBC-PCR assays [24]. However, some M. tuberculosis strains lacking the IS6110 gene were found in several clinical investigations [25, 26]. Thus, the mpb64 gene, which encodes the RD2 region of the MPB64 protein in the genome of M. tuberculosis [27], was introduced in our study, and it showed excellent specificity for MTBC strains in previous studies [9, 27]. There is some evidence for the absence of the mpb64 gene in some substrains of Mycobacterium bovis BCG [28]. However, there are no reports of the absence of both the IS6110 and mpb64 genes. In this study, we chose both the IS6110 gene and mpb64 gene as the target genes for MTBC to ensure its specificity. The MCDA-LFB technique targeting the *IS6110* and *mpb64* genes was successfully established and applied to the detection of MTBC. As Table 2 shows, 38 MTBC strains, 12 NTM strains and 10 other bacterial strains were used to test the specificity of the MCDA-LFB assay. The specificity of the MCDA-LFB assay was 100%, and there was no cross-reactivity with other pathogens (including NTM and nonmycobacterial strains).

In this research, the MCDA products were analyzed using the LFB, MG indicator and real-time turbidity methods. Comparing the above three validation methods, the LFBs were more convenient, reliable and visualized. They did not require special instruments or reagents [16, 17]. In addition, LFB can detect both the IS6110 and mpb64 genes and allow them to be visualized in a single test at the same time. According to the results, at 67 °C, both IS6110 and mpb64 gene amplification had higher turbidity and shorter reaction times at 30 min and 23 min (Fig. 2 and Table 1). The LFB assay indicated that the earliest test lines were observed for both TL1 and TL2 when the amplification lasted for 30 min at 67 °C, but the TL2 line was light red (Fig. 3). Hence, to ensure adequate amplification, the conditions of amplification temperature and time were optimized for MCDA-LFB at 67°C for 35 min in this study. In addition, based on the IS6110 and mpb64 genes, the MCDA assay LoD was 100 fg of genomic DNA per reaction (Fig. 4). We found that the sensitivity of the single MCDA assay (10 fg per reaction) was higher than that of the multiplex MCDA assay (100 fg per reaction) in the sensitivity experiment. This phenomenon may be caused by the different concentrations of the single set of primers in the singlex and/ or multiplex MCDA reaction system (e.g., the CP1 of IS6110 primers was 2.4μM in the singlex MCDA system and $1.2\,\mu\text{M}$ in the multiplex MCDA system) [29].

For sputum sample detection, Xpert MTB/RIF and MCDA-LFB had higher positive rates than the conventional culture and acid-fast staining assays (Table 3). While this is a possible reason for culture method that there were no living MTBC strains in the sputum samples, the Xpert MTB/RIF and/or MCDA-LFB could detect the nucleic acid of the target strains. In addition, the positive rates of the Xpert MTB/

(See figure on next page.)

Huang et al. BMC Microbiol (2021) 21:267 Page 6 of 11




Fig. 4 (See legend on previous page.)

Huang et al. BMC Microbiol (2021) 21:267 Page 7 of 11

Table 2 The details of the strains

Bacteria	Strain no. (source of strains)	No. of strains	MCDA- LFB result
MTBC			
M. tuberculosis	H37RvATCC27294	1	Р
M. tuberculosis	Isolated strains (GZCDC)	34	Р
M. bovis	ATCC19210	1	Р
M. africanum	ATCC25420	1	Р
Bacillus Calmette-Guerin	Vaccine strain (GZCDC-BCG)	1	Р
NTM			
M. aureus	l. aureus ATCC23366		N
M. flavum	ATCC43999	1	N
M. avium	ATCC25291	1	N
M. marinum	ATACC927	1	N
M. abscess	ATCC19977	1	N
M. chelonae	ATCC14472	1	N
M. gordon	ATCC14470	1	N
M. phlei	ATCC11758	1	N
M. nonchromogenic	ATCC19530	1	N
M. xenopi	ATCC19250	1	N
M. aichiense	ATCC27280	1	N
M.microti	ATCC19422	1	N
Other bacteria species			
Enterococcus faecalis	Isolate strains (GZCDC)	1	N
Salmonella	Isolate strains (GZCDC)	1	N
Klebsiella pneumoniae	Isolate strains (GZCDC)	1	N
Pseudomonas aeruginosa	Isolate strains (GZCDC)	1	N
Staphylococcus aureus	Isolate strains (GZCDC)	1	N
Escherichia coli	Isolate strains (GZCDC)	1	N
Bacillus cereus	Isolate strains (GZCDC)	1	N
Listeria monocytogenes	Isolate strains (GZCDC)	1	N
Streptococcus pneumoniae	Isolate strains (GZCDC)	1	N
Campylobacter jejuni	Isolate strains (GZCDC)	1	N

Abbreviations: MTBC Mycobacterium tuberculosis complex, NTM nontuberculous mycobacterium, ATCC American Type Culture Collection, GZCDC Guizhou Provincial Center for Disease Control and Prevention, P positive, N negative

RIF (76.47%, 39/51) and MCDA-LFB assays (80.39%, 41/51) for sputum specimens were similar (Table 3). Xpert MTB/RIF is an important tool for MTBC with

Table 3 Sputum samples and MTBC strain detection

Methods	Sputum samples (N = 51)	
	Positive	Negative
Acid-fast staining	26	25
Culture	35	16
Gene-Xpert	39	12
MCDA-LFB	41	10

Abbreviations: MCDA-LFB multiple cross displacement amplification with lateral flow biosensor

rifampin (RIF) resistance [30], and the detection time is approximately 2.5 h. Usually, Xpert MTB/RIF has a four-channel detection system, which means that this machine can test 4 samples at a time and is not suitable for the large-scale detection of MTB samples. Moreover, Xpert MTB/RIF needs to use the matching kit, which represents a high cost. However, the MCDA-LFB technique was a better choice for detecting MTBC with a large specimen size. Importantly, it is low cost (approximately 7.5 USD per test). Compared with the AFS method, combined microscopy detected MTB of sputum specimens, an approach that lacks sensitivity [31, 32] and takes approximately 2 h, including fixing, dyeing, drying naturally and microscopy. Microscopy

Huang et al. BMC Microbiol (2021) 21:267 Page 8 of 11

took more time as the samples increased. However, the MCDA-LFB technique can solve this problem well. The whole detection process was completed within 70 min, including genomic template preparation (approximately 30 min), MCDA reaction (approximately 35 min) and LFB verification (approximately 1–2 min), and the maximum detection quantity was 96 samples with a conventional PCR apparatus.

Conclusions

In this study, a reliable, rapid, visual, inexpensive and simple MCDA-LFB assay based on the *IS6110* and *mpb64* genes was successfully established for the detection of MTBC. This technique could identify the *M. tuberculosis* complex with high specificity, sensitivity, and rapidity, and it was a visual and simple assay of the MCDA products. Thus, the MTBC-MCDA-LFB method could be regarded as a useful technique for rapid and reliable identification of the *M. tuberculosis* complex in clinical samples. Moreover, the technique did not require complicated instruments and expensive reagents and could be used widely, especially in resource-limited areas of developing countries with a high TB burden.

Methods

Reagents and apparatus

DNA isothermal amplification kits, colorimetric indicators (malachite green, MG), and biotin-14-dCTP were purchased from Bei-Jing HaiTaiZhengYuan. Co., Ltd. (Beijing, China). DNA extraction kits (QIAamp DNA minikits; Qiagen, Hilden, Germany) were purchased from Qiagen (Beijing, China). The nanoparticle-based lateral flow biosensor (Disposable Lateral Flow Biosensor, LFB) was manufactured by Tian-Jin HuiDeXin Biotech. Co., Ltd. (Tianjin, China).

Design of MCDA primers

Two specific MTC-MCDA target genes, *IS6110* (Gen-Bank, Sequence ID: CP053903.1) and *mpb64* (Gen-Bank, Sequence ID: CP053903.1), were chosen. The 10 specific primers were designed by primer software PRIMER PREMIER 5.0 and Primer Explorer V4 [15]. In addition, FITC (fluorescein isothiocyanate) was labeled at the 5' end of the *IS6110*-C1* primer, and digoxigenin (Dig) was labeled at the 5' end of the *mpb64*-C1* primer. The details, including primer design, sequences and modifications, and locations in the expression sites of the *IS6110* and *mpb64* genes are listed in Table 4 and Fig. 5. All of the primer sequences were synthesized

Table 4 The details of primers for the *IS6110* and *mpb64* genes

Genes	Primers name ^a	Sequences and modifications	Length
IS6110	IS6110-F1	5'-GGATGGTCGCAGAGATCC-3'	18 nt
	IS6110-F2	5'-ATCGCGTTCGCCCTT-3'	15 nt
	IS6110-CP1	5'-CGCGCAGCCAACACCAAGTAGCAGCACGATTCGGAGTG-3'	38 mer
	IS6110-CP2	5'-CCGGGACCACGACCGAAGACGCAATTCGGCGTTGTC-3'	36 mer
	IS6110-C1*	5'-FITC-CGCGCAGCCAACACCAAGTAG-3'	21 nt
	IS6110-C2	5'-CCGGGACCACGACCGAAGA-3'	19 nt
	IS6110-D1	5'-ACCTCACTGATCGCTG-3'	16 nt
	IS6110-D2	5'-ATCCGCTGAGCTGAAGC-3'	17 nt
	IS6110-R1	5'-ACTTACGCACCGTCTC-3'	16 nt
	IS6110-R2	5'-CAGGCGCAGGTCGATG-3'	16 nt
mpb64	mpb64-F1	5'-CCCCGGGTTGAAGAAGA-3'	17 nt
	mpb64-F2	5'-GCTCAAGGTCTACCAGAAC-3'	19 nt
	mpb64-CP1	5'-ACAGGTATCGATAGCGCCGAATGCCCCGTCGTTCGTGACT-3'	40 mer
	mpb64-CP2	5'-TGCCACAGCGTGTCATAGGTACGACCACGTACAAGGC-3'	37 mer
	mpb64-C1*	5'-Dig-ACAGGTATCGATAGCGCCGAATG-3'	23 nt
	mpb64-C2	5'-TGCCACAGCGTGTCATAGGT-3'	20 nt
	mpb64-D1	5'-CGGTGAATTATCAGAACTTC-3'	20 nt
	mpb64-D2	5'-GCCTGGTCCCAATCGAA-3'	17 nt
	mpb64-R1	5'-GTGAACTGAGCAAGCAGA-3'	18 nt
	mpb64-R2	5'-ACAATGGGGAAGACGACT-3'	18 nt

Abbreviations: FITC fluorescein isothiocyanate, Dig digoxigenin, mer monomeric unit, nt nucleotide

a IS6110-C1*, 5'-labeled with FITC when used in the MCDA-LFB assay; mpb64-C1*, 5'-labeled with Dig when used in the MCDA-LFB assay

Huang et al. BMC Microbiol (2021) 21:267 Page 9 of 11




Fig. 5 Sequences, modifications and locations in the expression sites of the *IS6110* and *mpb64* genes. **A** Location and nucleotide sequence of the *M. tuberculosis IS6110* gene used to design the MCDA primers. **B** Location and nucleotide sequence of the *M. tuberculosis mpb64* gene used to design the MCDA primers. All sequences of the primer sites are underlined. Right and left arrows indicate sense and complementary sequences that are used

and purified through Tian-Yi Biotech (Beijing, China) at HPLC purification grade.

Bacterial strains and genomic DNA preparation

There were 60 strains, including 38 MTBC strains, 12 nontuberculous mycobacterium (NTM) strains and 10 other bacteria, and the details are shown in Table 2. The bacterial strains were stored in 10% (w/v) glycerol broth at $-80\,^{\circ}\text{C}$ and then revived and cultured. According to the instruction book, the genomic templates were extracted from the cultured strains by the QIAamp DNA Mini Kit (Qiagen, Germantown, MD, USA). Sputum samples were treated with 4% NaOH solution, and genomic templates were extracted by the kit. Then, the genomic templates were tested by ultraviolet spectrophotometry at A260/280 and stored at $-20\,^{\circ}\text{C}$ before the templates were used.

MCDA reaction and detection

The suitability of two target genes for the MCDA primers was examined through MCDA reaction of the single IS6110 gene for the mpb64 gene. Then, both target genes were detected at the same time. Genomic templates of the standard strain M. tuberculosis (H37Rv, ATCC 27294) were used for the MCDA assay according to the DNA isothermal amplification kits. Briefly, the MCDA reaction mixtures contained 12.5 µl of $2 \times$ Buffer, 1 µl of Bst 2.0 DNA polymerase (10 U), 1 µl of biotin-14-dCTP, 1 µl of malachite green, 2 µl of genomic templates, 5.3 µl of DW and 2.2 µl of amplification primers or mixed primers containing 0.4 µM each of displacement primers F1 and F2, 0.4 µM each of amplification primers C1* and C2, 1.2 µM each of amplification primers R1, R2, D1 and D2, and 2.4 µM each of cross primers CP1 and CP2 in a final volume of 25 µl. Then, Huang et al. BMC Microbiol (2021) 21:267 Page 10 of 11

they were reacted for $35\,\mathrm{min}$ at $67\,^\circ\mathrm{C}$. Three detection methods, including a real-time turbidimeter (LA-320C), disposable lateral flow biosensor (LFB) and colorimetric indicator (MG), were applied to detect the MCDA amplification products.

The threshold value (turbidity) was 0.1, and a turbidity of > 0.1 was regarded as positive amplification for the MCDA-LFB assay by using a real-time turbidimeter [16]. The LFB consisted of the sample pad, conjugate pad, NC membrane, and absorbent pad, which were laminated onto a plastic adhesive backing card. Then, the anti-FITC Ab, anti-Dig Ab and biotin-BSA were sprayed onto the NC membrane for TL1, TL2 and CL, respectively [19]. The detection method of LFB involves depositing an aliquot (1-2 µl) of MCDA amplification products on the sample pad of LFB and then depositing an aliquot of running buffer (150–200 µl) on the sample pad of LFB. The amplicons (crimson red) were captured by specific antibodies (namely, anti-FITC Ab, anti-Dig Ab, and biotin-BSA) when flowing through the NC membrane. The red lines (i.e., CL, TL1 and/or TL2) were observed in the positive amplification, while only CL was red in the blank and negative controls. In addition, the color of positive amplification changed from blue to light blue, while the blank and negative controls were colorless by using MG indicators. In this study, DW was used as the template for the blank control, and the genomic DNA of M. avium was used as the template in the negative control.

Optimizing the reaction temperature and time of the MCDA assay

Both the *IS6110* gene and *mpb64* gene were detected at different temperatures by a real-time turbidimeter (LA-320C). We examined the effect of different temperatures, from 63 to 70°C, with 1°C intervals for MCDA amplification. The *M. avium* genomic templates were used as a negative control. Subsequently, the two target genes were confirmed at different amplification reaction times, 20, 30, 40 and 50 min, detected by LFBs and MG.

Analytical sensitivity and specificity of MTBC-MCDA-LFB assay

The genomic DNA templates of MTB (H37Rv, ATCC 27294) were serially diluted (100 ng-1 fg per microliter) for sensitivity analysis by MCDA-LFB detection for single genes and double genes. DW was the template for the blank control. Thirty-eight MTBC strains, 12 NTM strains and 10 other bacteria (Table 2) were used for the specificity assay. The genomic DNA from the strains was amplified by MCDA reactions and then

assayed through LFB. The experiments were repeated at least two times.

Practicability of the MTBC-MCDA-LFB assay for sputum samples

Fifty-one sputum samples provided by the Pulmonary Hospital of Guiyang (China) were detected by AFS, conventional culture, Xpert MTB/RIF and MTBC-MCDA-LFB methods. The AFS method was performed in our study according to previous publications [32]. Meanwhile, Xpert MTB/RIF and conventional culture methods were implemented according to the manufacturer's instructions. The MTBC-MCDA assay was carried out as described above, and its amplicons were validated by LFBs.

Abbreviations

MTBC: Mycobacterium tuberculosis complex; TB: Tuberculosis; MCDA: Multiple cross displacement amplification; LFB: Nanoparticle-based lateral flow biosensor; AFS: Acid-fast staining; WHO: World Health Organization; FITC: Fluorescein isothiocyanate; Dig: Digoxigenin. mer, monomeric unit; nt: Nucleotide; ATCC: American Type Culture Collection; LoD: Limit of detection; MG: Malachite green; GZCDC: Guizhou Provincial Center for Disease Control and Prevention; NTM: Nontuberculous mycobacterium; P: Positive; N: Negative; TL1: Test Line 1; TL2: Test Line 2; CL: Control line; DW: Double-distilled water; BC: Blank control; PC: Positive control; NC: Negative control.

Acknowledgments

We thank all of the medical workers in both the Public Health School, Guizhou Medical University and The Second Affiliated Hospital, Guizhou University of Traditional Chinese Medicine, for their cooperation in this study.

Authors' contributions

JH and SL conceived and designed the study. JH, ZX and SL participated in primer design. JH, ZX, XY, and XC contributed to all of the laboratory work. WX and YC contributed to the data collection. WZ, WC and HC performed the statistical analysis. JH and SL wrote the initial draft of the manuscript. JH and SL revised the manuscript. All authors read and approved the final manuscript.

Fundina

This study was funded by grants from the Science and Technology Department of Guizhou Province (No. [2019]2822, [2019]1186), [2018]-5606 and [2018]-5767).

Availability of data and materials

The datasets used and/or analyzed during the current study are available from the corresponding author upon reasonable request. The raw sequence data reported in this paper came from GenBank, Sequence ID: CP053903.1, *M. tuberculosis* strain H37Rv_IC1chromosome (https://www.ncbi.nlm.nih.gov/nuccore/CP053903.1?report=fasta).

Declarations

Ethics approval and consent to participate

The study was approved by the Human Ethics Committee of the Guizhou Provincial Center for Disease Control and Prevention and complied with the Declaration of Helsinki. All data/isolates were analyzed anonymously.

Consent for publication

Not applicable.

Competing interests

All of the authors declare that there are no competing interests in this article.

Huang et al. BMC Microbiol (2021) 21:267 Page 11 of 11

Author details

¹Laboratory of Infectious Disease of Experimental Center, Guizhou Provincial Center for Disease Control and Prevention, 73 Bageyan Road, Guiyang 550004, Guizhou, People's Republic of China. ²Public Health School, Guizhou Medical University, Guiyang 550025, Guizhou, People's Republic of China. ³The Second Affiliated Hospital, Guizhou University of Traditional Chinese Medicine, Guiyang, Guizhou 550003, People's Republic of China.

Received: 7 December 2020 Accepted: 24 September 2021 Published online: 04 October 2021

References

- Fieweger RA, Wilburn KM. VanderVen BC: comparing the metabolic capabilities of Bacteria in the *Mycobacterium tuberculosis* complex. Microorganisms. 2019;7(6).
- Cadmus SI, Jenkins AO, Godfroid J, Osinusi K, Adewole IF, Murphy RL, et al. Mycobacterium tuberculosis and Mycobacterium africanum in stools from children attending an immunization clinic in Ibadan, Nigeria. Int J Infect Dis. 2009:13(6):740–4.
- Rodriguez-Campos S, Smith NH, Boniotti MB, Aranaz A. Overview and phylogeny of Mycobacterium tuberculosis complex organisms: implications for diagnostics and legislation of bovine tuberculosis. Res Vet Sci. 2014;97(Suppl):S5–S19.
- Tientcheu LD, Koch A, Ndengane M, Andoseh G, Kampmann B, Wilkinson RJ. Immunological consequences of strain variation within the *mycobacterium tuberculosis* complex. Eur J Immunol. 2017;47(3):432–45.
- Organization WHO: Global tuberculosis report 2019. Online at https:// www.who.int/teams/global-tuberculosis-programme/tb-reports/ global-report-2019.
- Mwaba P, McNerney R, Grobusch MP, O'Grady J, Bates M, Kapata N, et al. Achieving STOP TB partnership goals: perspectives on development of new diagnostics, drugs and vaccines for tuberculosis. Tropical Med Int Health. 2011;16(7):819–27.
- Zumla A, Atun R, Maeurer M, Mwaba P, Ma Z, O'Grady J, et al. Viewpoint: scientific dogmas, paradoxes and mysteries of latent *Mycobacterium* tuberculosis infection. Trop Med Int Health. 2011;16(1):79–83.
- Schito M, Migliori GB, Fletcher HA, McNerney R, Centis R, D'Ambrosio L, et al. Perspectives on advances in tuberculosis diagnostics, drugs, and vaccines. Clin Infect Dis. 2015;61(suppl 3):S102–18.
- Lekhak SP, Sharma L, Rajbhandari R, Rajbhandari P, Shrestha R, Pant B. Evaluation of multiplex PCR using MPB64 and IS6110 primers for rapid diagnosis of tuberculous meningitis. Tuberculosis (Edinb). 2016;100:1–4.
- Ji L, Jiang Y, Li G, Zhao X, Wan K. A real-time PCR assay based on a specific mutation of *PstS1* gene for detection of *M. bovis* strains. Biologicals. 2020;64:23–7.
- Sinha P, Prakash P, Patne SC, Anupurba S, Gupta S, Srivastava GN. Performance of nested multiplex PCR assay targeting MTP40 and IS6110 gene sequences for the diagnosis of tubercular lymphadenitis. J Microbiol. 2017;55(1):63–7.
- Chakravorty S, Simmons AM, Rowneki M, Parmar H, Cao Y, Ryan J, et al. The new Xpert MTB/RIF ultra: improving detection of *Mycobacterium tuberculosis* and resistance to rifampin in an assay suitable for point-of-care testing. MBio. 2017;8(4):e00812–17.
- Joon D, Nimesh M, Varma-Basil M, Saluja D. Evaluation of improved IS6110 LAMP assay for diagnosis of pulmonary and extra pulmonary tuberculosis. J Microbiol Methods. 2017;139:87–91.
- Fang R, Li X, Hu L, You Q, Li J, Wu J, et al. Cross-priming amplification for rapid detection of mycobacterium tuberculosis in sputum specimens. J Clin Microbiol. 2009;47(3):845–7.
- Wang Y, Wang Y, Ma A-J, Li D-X, Luo L-J, Liu D-X, et al. Rapid and sensitive isothermal detection of nucleic-acid sequence by multiple cross displacement amplification. Sci Rep. 2015;5:11902.
- Li S, Liu C, Liu Y, Ma Q, Wang Y. Development of a multiple cross displacement amplification combined with nanoparticles-based biosensor assay to detect Neisseria meningitidis. Infect Drug Resist. 2019;12:2077–87.

- Li S, Liu Y, Wang Y, Wang M, Liu C. Rapid detection of *Brucella spp.* and elimination of carryover using multiple cross displacement amplification coupled with nanoparticles-based lateral flow biosensor. Front Cell Infect Microbiol. 2019;9:78.
- Wang Y, Li H, Wang Y, Li H, Luo L, Xu J, et al. Development of multiple cross displacement amplification label-based gold nanoparticles lateral flow biosensor for detection of *listeria monocytogenes*. Int J Nanomedicine. 2017;12:473–86.
- Wang Y, Wang Y, Zhang L, Xu J, Ye C. Visual and multiplex detection of nucleic acid sequence by multiple cross displacement amplification coupled with gold nanoparticle-based lateral flow biosensor. Sensors Actuators B Chem. 2017;241:1283–93.
- Wang Y, Hui L, Li D, Li K, Wang Y, Xu J, et al. Multiple cross displacement amplification combined with gold nanoparticle-based lateral flow biosensor for detection of vibrio parahaemolyticus. Front Microbiol. 2016;7(770):2047
- 21. Waterfield T, Fairley D, Lynn F, Blackwood B, Shields MD. A protocol for a systematic review of the diagnostic accuracy of loop-mediated-isothermal AMPlification (LAMP) in diagnosis of invasive meningococcal disease in children. Syst Rev. 2018;7(1):86.
- Zhang W, Chen C, Cui J, Bai W, Zhou J. Application of loop-mediated isothermal amplification (LAMP) assay for the rapid diagnosis of pathogenic bacteria in clinical sputum specimens of acute exacerbation of COPD (AECOPD). Int J Clin Exp Med. 2015;8(5):7881–9.
- 23. Yang Q, Domesle KJ, Wang F, Ge B. Rapid detection of *Salmonella* in food and feed by coupling loop-mediated isothermal amplification with bioluminescent assay in real-time. BMC Microbiol. 2016;16(1):112.
- Thierry D, Eisenach KD, Crawford JT, Bates JH, Gicquel B, Guesdon JL. IS6110, an IS-like element of Mycobacterium tuberculosis complex. Nucleic Acids Res. 1990;18(1).
- Sharma K, Sharma M, Batra N, Sharma A, Dhillon MS. Diagnostic potential of multi-targeted LAMP (loop-mediated isothermal amplification) for osteoarticular tuberculosis. J Orthop Res. 2017;35(2):361–5.
- Víctor Vinuesa RB, Briones ML, Clari MÁ, Cresencio V, Giménez E, Carmen Muñoz RO, et al. Performance of a highly sensitive mycobacterium tuberculosis complex real-time PCR assay for diagnosis of pulmonary tuberculosis in a low-prevalence setting: a prospective intervention study. J Clin Microbiol. 2018;56(5):e00116–8.
- Tamada Y, Kanda S, Yoshidome A, Hayashi I, Miyake M, Nishiyama T. Diagnosis of active tuberculosis using MPB64, a specific antigen of Mycobacterium bovis. Microbiol Immunol. 2012;56(11):740–7.
- Huayi LI, Ulstrup JC, Jonassen TO, Melby K, Nagai S, Harboe M. Evidence for absence of the MPB64 gene in some substrains of Mycobacterium bovis BCG. Infect Immun. 1993;61(5):1730–4.
- Yang X, Huang J, Chen X, Xiao Z, Wang X, Chen Y, et al. Rapid and visual differentiation of mycobacterium tuberculosis from the mycobacterium tuberculosis complex using multiplex loop-mediated isothermal amplification coupled with a nanoparticle-based lateral flow biosensor. Front Microbiol. 2021;12:708658.
- Boehme MD CC, Pamela Nabeta MD, Hillemann D, Nicol MP, Shubhada Shenai PD, Krapp F, et al. Rapid molecular detection of tuberculosis and rifampin resistance. N Engl J Med. 2010;363(11):1005–15.
- Raj A, Singh N, Gupta KB, Chaudhary D, Yadav A, Chaudhary A, et al. Comparative evaluation of several gene targets for designing a multiplex-PCR for an early diagnosis of Extrapulmonary tuberculosis. Yonsei Med J. 2016:57(1):88
- Agrawal M, Bajaj A, Bhatia V, Dutt S. Comparative study of GeneXpert with ZN stain and culture in samples of suspected pulmonary tuberculosis. J Clin Diagn Res. 2016;10(5):DC09–12.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.